



APPLICATION INSTRUCTIONS

(Strap may contain latex always ask patient if they are allergic)

| | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| <p>Step 1: Carefully unwrap the Shower Boot. Slide the Shower Boot over leg to be protected. Slightly twist Shower Boot at top to bring the strap from behind the leg to the front of the leg.</p> | |
| <p>Step 2. Fold excess plastic and begin to wrap the strap. (<i>Make sure no plastic is showing above the strap, as the strap is the water barrier.</i>)</p> | |
| <p>Step 3: As you wrap around the leg, PULL strap tightly to keep the strap snug and ensure closure of the Shower Boot to the leg. Lightly squeeze strap around leg when finished. (<i>CAUTION: Plastic may be slippery when wet</i>)</p> | |
| <p>Step 4: After showering, dry the Shower Boot and surrounding area. Carefully remove the Shower Boot by loosening the strap. Dampness due to perspiration is normal, just pat dry.</p> | |