




APPLICATION INSTRUCTIONS

(Strap may contain latex always ask patient if they are allergic)

<p>Step 1: APPLY AND GUIDE: Carefully unwrap the Shower Glove. Slide the glove over the hand and arm to be protected. Have patient fan fingers. Push down glove in between each finger so the glove fits securely on the hand and fingers.</p>	
<p>Step 2: TWIST AND FOLD: Lay strap over top and slightly TWIST glove to the outside of arm. <i>(For a smaller arm, at wrist, continue to twist up and out tightening the glove securely around the full length of arm.)</i> Fold excess plastic to the outside at top of glove before beginning to wrap.</p>	
<p>Step 3 : WRAP AND SQUEEZE <i>(The strap is the water barrier, make sure no plastic is showing above the strap.)</i> As you wrap, PULL strap tightly to keep the strap snug. When finished wrapping, lightly squeeze strap to ensure closure of the glove to the arm.</p>	
<p>Step 4: After showering, dry the Shower Glove and surrounding area. Carefully remove the glove by loosening the strap. Dampness due to perspiration is normal, just pat dry.</p>	



www.ShowerCovers.com